


Single Nucleus Sequencing Enabled by the CellRaft Technology



CELL
Microsystems

Nick Trotta, PhD
Cell Microsystems, Inc.



ntrotta@cellmicrosystems.com

CellRaft Technology


- Straightforward, “track and trace” means of imaging, sorting and isolating single cells or nuclei
- Employs a microwell array which replicates standard culture conditions – improves cell viability and reduces appearance of spurious phenotypes
- Multiple products available; each tailored to specific experimental and throughput needs

CellRaft Technology - Components

- **CytoSort Array**
 - Microwell array
 - Releasable CellRaft in each microwell
 - 100 or 200 micron CellRafts
 - Elastomer array
 - Polystyrene CellRaft
- **Release/Transfer Devices**
 - Needle for release of CellRaft from array
 - Magnetic wand for CellRaft recovery and transfer to plate or tube
 - Manual and automated versions



CellRaft Technology - Operation



CellRaft Systems


■ CellRaft System – Manual

- Fits on virtually any inverted microscope
- Very flexible and cost-effective
- Compatible with any type/size of CytoSort Array
- Approximately 50 cells in 2–4 hours
- Manual collection into any vessel
- Comes with 2 CytoSort Arrays



■ AIR™ System – Automated

- Fully integrated microscope and release/transfer hardware
- 3-channel fluorescence
- “Real-time” or “Cytometric” sorting capabilities
- 96 cells in approximately 1.5 hr;
- 30 min for each subsequent 96 cells
- Collect into
 - 96-well plates (~150–200 μ L volume)
 - PCR strip tubes (minimum 2.5 μ L volume).



Contact Us

Contact us for info and materials

- info@cellmicrosystems.com
- ntrotta@cellmicrosystems.com



CellRaft AIR™ System

-  **Fully automated**
Imaging, sorting and isolation of single cells with *no minimum sample size requirement*
-  **Application flexibility**
Compatible with genomics and cloning workflows
-  **Standard cell culture**
CytoSort Arrays keep cells happy and viable, enabling cloning and reducing artifacts
-  **Imaging-based sorting**
Detailed image analysis with 3-channel fluorescence
-  **1 hour for 96 cells**
Cells sorted into standard multi-well plates and tubes

info@cellmicrosystems.com



RaftNote

Applications and Protocols for Use with CellRaft Technology



Methods for Preparing Cryopreserved Neural Tissue Samples for Single Nucleus Sequencing (SNS) Using the CellRaft™ Technology

Mike McConnell, PhD¹, Jacquelyn DuVall, PhD², Nick Trotta, PhD².

1. University of Virginia; Department of Biochemistry and Molecular Genetics; 2. Cell Microsystems, Inc.; Research Triangle Park, NC

Abstract

Single cell genomic analysis has already impacted a range of life science fields, including tumor clonality, neuronal mosaicism, and drug resistance mechanisms. These methods rely on access to relatively fresh samples, ensuring that a maximum amount of high quality nucleic acid can be extracted from a given cell. Conversely, translational research depends on the availability of human samples which are difficult to acquire in large numbers as fresh samples. As an alternative to single cell sequencing, single nucleus sequencing (SNS) has enabled genomic analysis using the plentiful cryopreserved samples found in various biobanks and other facilities. The CellRaft Technology uniquely enables the sorting and isolation of single nuclei for downstream analysis. Here we provide optimal methods for imaging nuclei on CytoSort Arrays, and downstream molecular biology protocols for the amplification and preparation of nucleic acid prior to sequencing.

Translational Single Cell Analysis using SNS

Translational research relies on patient samples which are often stored in biobank facilities. Nucleic acid extraction from preserved samples often results in compromised integrity rendering single cell analysis challenging. To address this shortcoming, many investigators are turning to single nucleus isolation prior to molecular analysis. Here, we provide details on methods for isolating nuclei from cryopreserved neural tissue samples for

subsequent genomic analysis. The methods for SNS described here using the CellRaft Technology were developed by the laboratory of Mike McConnell, PhD of the University of Virginia. [6] The workflow steps are shown in Figure 1. Once nuclei from specific cell types have been isolated (in this case neurons), the nuclei are lysed and the DNA is amplified and prepared for sequencing. The CellRaft Technology allows nuclei to be isolated and sorted for cell type-specific markers prior to downstream molecular analysis.





Figure 1: Schematic of the SNS whole genome sequencing workflow using CellRaft Technology.





Single Nucleus Sequencing Enabled by CellRaft Technology

MIKE McCONNELL
@mikemc43

The DEPARTMENT of BIOCHEMISTRY *and* MOLECULAR GENETICS




UNIVERSITY
of **VIRGINIA**


SCHOOL OF MEDICINE



How do human genomes encode human brains?



Human induced Pluripotent Stem Cells



Single Cell Genomics


METHOD OF THE YEAR



"We are just beginning to understand the molecular diversity of cells in the brain," says Thomas Insel, director of the US National Institute of Mental Health. "Single-cell methods will be critical, not only for defining the taxonomy of neurons and glia but for revealing the effects of experience or development on profiles of expression within a brain region."

NATURE METHODS | VOL.11 NO.1 | JANUARY 2014 | 13

Mosaic Copy Number Variation in Human Neurons



McConnell, et al. (2013) *Science*

McConnell Lab

Human Neurobiology • Single Cell Genomics


www.mcconnell-lab.org

The DEPARTMENT of BIOCHEMISTRY and MOLECULAR GENETICS



UNIVERSITY
of VIRGINIA

SCHOOL OF MEDICINE



**nature
medicine**

Chromosome instability is common in human cleavage-stage embryos


Evelyne Vanneste^{1,2,9}, Thierry Voet^{1,9}, Cédric Le Caignec^{1,3,4}, Michèle Ampe⁵, Peter Konings⁶, Cindy Melotte¹, Sophie Debrock², Mustapha Amyere⁷, Miikka Vikkula⁷, Frans Schuit⁸, Jean-Pierre Fryns¹, Geert Verbeke⁵, Thomas D’Hooghe², Yves Moreau⁶ & Joris R Vermeesch¹

LETTER

doi:10.1038/nature09807

Tumour evolution inferred by single-cell sequencing

Nicholas Navin^{1,2}, Jude Kendall¹, Jennifer Troge¹, Peter Andrews¹, Linda Rodgers¹, Jeanne McIndoo¹, Kerry Cook¹, Asya Stepansky¹, Dan Levy¹, Diane Esposito¹, Lakshmi Muthuswamy³, Alex Krasnitz¹, W. Richard McCombie¹, James Hicks¹ & Michael Wigler¹



ARTICLE

Received 15 Sep 2015 | Accepted 12 Feb 2016 | Published 19 Apr 2016 | Updated 14 Jun 2016

DOI: [10.1038/ncomms11022](https://doi.org/10.1038/ncomms11022)

OPEN

Nuclear RNA-seq of single neurons reveals molecular signatures of activation

Benjamin Lacar^{1,*}, Sara B. Linker^{1,*}, Baptiste N. Jaeger^{1,*}, Suguna Rani Krishnaswami², Jerika J. Barron¹, Martijn J.E. Kelder¹, Sarah L. Parylak¹, Apuã C.M. Paquola¹, Pratap Venepally², Mark Novotny¹, Carolyn O'Connor¹, Conor Fitzpatrick¹, Jennifer A. Erwin¹, Jonathan Y. Hsu¹, David Husband¹, Michael J. McConnell³, Roger Lasken² & Fred H. Gage¹

McConnell Lab


Human Neurobiology • Single Cell Genomics

www.mcconnell-lab.org

The DEPARTMENT of BIOCHEMISTRY and MOLECULAR GENETICS




SCHOOL OF MEDICINE




Saline versus PTZ

Intact cells lose transcriptional signature



Nuclei retain transcriptional signature

a






Read Depth


PE 50bp reads mapped to ~500kb non-overlapping bins of mapable sequence

Segmentation algorithm calls regions of copy number gain / loss based on MAD scores




(Note: not to scale)

Whole Genome Amplification (WGA) followed by single cell sequencing




Integer-like gain and
loss in trisomic male
fibroblasts

A Chr2 Duplication in a FCTX Neuron





A Chr1 Deletion in a FCTX Neuron



McConnell Lab

Human Neurobiology • Single Cell Genomics

www.mcconnell-lab.org


The DEPARTMENT of BIOCHEMISTRY and MOLECULAR GENETICS




UNIVERSITY
of VIRGINIA

SCHOOL OF MEDICINE


PicoPLEX
(Rubicon Genomics)




MALBAC



MDA



Low median absolute deviation (MAD) indicates high confidence CNV detection.
Ginkgo (CSHL)




McConnell Lab

Human Neurobiology • Single Cell Genomics

www.mcconnell-lab.org

The DEPARTMENT of BIOCHEMISTRY and MOLECULAR GENETICS



McConnell Lab

Human Neurobiology • Single Cell Genomics


www.mcconnell-lab.org

The DEPARTMENT of BIOCHEMISTRY and MOLECULAR GENETICS



UNIVERSITY
of VIRGINIA

SCHOOL OF MEDICINE



McConnell Lab

Human Neurobiology • Single Cell Genomics



www.mcconnell-lab.org

The DEPARTMENT of BIOCHEMISTRY and MOLECULAR GENETICS



UNIVERSITY
of VIRGINIA

SCHOOL OF MEDICINE



McConnell Lab

Human Neurobiology • Single Cell Genomics


www.mcconnell-lab.org

The DEPARTMENT of BIOCHEMISTRY and MOLECULAR GENETICS



UNIVERSITY
of VIRGINIA

SCHOOL OF MEDICINE




McConnell Lab

Human Neurobiology • Single Cell Genomics

www.mcconnell-lab.org

The DEPARTMENT of BIOCHEMISTRY and MOLECULAR GENETICS



McConnell Lab

Human Neurobiology • Single Cell Genomics

www.mcconnell-lab.org


The DEPARTMENT of BIOCHEMISTRY and MOLECULAR GENETICS




UNIVERSITY
of VIRGINIA

SCHOOL OF MEDICINE

MALBAC
(Standard Buffer)



MALBAC
(Burbulis Buffer)



Advantages of CellRaft Approach:

- Avoid using core facility
 - small number of cells required
 - rapid turnaround to next experiment
 - schedule as needed
- Positive identification of specific single cells
- Preserve scarce tissue
 - >50,000 cells needed to establish FACS gates

McConnell Lab

Human Neurobiology • Single Cell Genomics

www.mcconnell-lab.org

The DEPARTMENT of BIOCHEMISTRY and MOLECULAR GENETICS



Single Cell Analysis

June 29 - July 12, 2018

Application Deadline: March 31, 2018

Instructors:

David Chenoweth, University of Pennsylvania
Michael McConnell, University of Virginia School of Medicine
Gene Yeo, University of California, San Diego

Co-Instructor:

See the [roll of honor](#), who's taken the course in the past

The goal of this two-week course is to familiarize students with cutting-edge technologies for characterization of single cells. Modules of the course will be taught by scientists with expertise in distinct areas of single cell analysis. Topics to be covered include quantitative single cell analysis by RNAseq, genomic DNA analysis, proteomics, and metabolomics. Multiple nucleic acid amplification methodologies including droplet-based RNAseq, MALBAC and MDA will be employed. In addition, students will be instructed in basic bioinformatic analysis of next generation sequencing data.



Topics:

- Single cell genome, transcriptome, and proteome measurement
- Introductory next generation sequencing data analysis
- Photoactivatable single cell probes
- Single cell mass spectrometry / Soft X-ray tomography

Speakers and Module Leaders in 2017 include:

Nancy Albritton, University of North Carolina
Lacramioara Bintu, Stanford University
Jim Eberwine, University of Pennsylvania
Amy Herr, University of California Berkeley
Carolyn Larabell, University of California, San Francisco
Elena Romanova, University of Illinois
Stas Rubakhin, University of Illinois
Rickard Sandberg, Karolinska Institute, Sweden
Carsten Schultz, Oregon Health & Science University
Peter Sims, Columbia University
Nick Trotta, Cell Microsystems
Xiaowei Zhang, Harvard University

This course is supported with funds provided by: [National Institute of General Medical Sciences](#), [Howard Hughes Medical Institute](#), and [Helmsley Charitable Trust](#).



Acknowledgements

www.McConnell-Lab.org @mikemc43



Ian Burbulis, Ph.D.
Maggie Wierman, Ph.D.
Lise Harbom
Will Chronister
Nadine Michel
Matt Wolpert
Mark Haakenson

Collaborators
Rusty Gage (Salk)
Dan Weinberger (LIBD)
Kristen Brennand (Mt. Sinai)
Mark Beenhakker (UVa)
Aakrosh Ratan (UVa)
Stefan Bekiranov (UVa)

Funding: NIMH U01 MH106882